




# Application of the Original von Rosen splint



*The hips must be flexed more than 90 degrees and abducted 60-70 degrees.  
Some hip motion must be possible in the splint and forced abduction be avoided.*

Manufacturer of the Original von Rosen splint

ISAKSSONS GUMMIFABRIK AB

Ängalagsvägen 226,

SE-269 95 Båstad

Sweden

Phone: +46 431 36 31 60 Fax: +46 431 36 31 82

E-mail: [isakssons@bjarenet.com](mailto:isakssons@bjarenet.com)

[www.originalvonrosen.com](http://www.originalvonrosen.com)

Copyright©Medi-Estate 2013